

Hydromorphone Hydrochloride Oral Solution

DEFINITION

Hydromorphone Hydrochloride Oral Solution contains NLT 90.0% and NMT 110.0% of the labeled amount of hydromorphone hydrochloride ($C_{17}H_{19}NO_3 \cdot HCl$). It may contain suitable preservatives.

IDENTIFICATION

- A.** The retention time of the major peak of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in the *Assay*.

ASSAY

PROCEDURE

Diluent: Phosphoric acid and water (1:1000)

Solution A: 1.0 mg/mL of sodium 1-heptanesulfonate monohydrate in methanol and water (1:9). To each liter of this solution add 1.0 mL of triethylamine, and adjust with phosphoric acid to a pH of 2.5 ± 0.1 .

Solution B: 1.0 mg/mL of sodium 1-heptanesulfonate monohydrate in methanol and water (1:1). To each liter of this solution add 1.0 mL of triethylamine, and adjust with phosphoric acid to a pH of 2.5 ± 0.1 .

Mobile phase: See *Table 1*.

Table 1

Time (min)	Solution A (%)	Solution B (%)
0	85	15
24	5	95
25	85	15
30	85	15

The *Standard solution* and *Sample solution* should be kept in a cool place, protected from light.

Standard solution: 0.08 mg/mL of USP Hydromorphone Hydrochloride RS in *Diluent*

Sample solution: Nominally 0.08 mg/mL of hydromorphone hydrochloride obtained by diluting a suitable volume of Oral Solution in *Diluent*

Chromatographic system

(See *Chromatography* (621), *System Suitability*.)

Mode: LC

Detector: UV 220 nm

Column: 4.6-mm \times 5-cm; 3.5- μ m packing L1

Column temperature: 45°

Flow rate: 1 mL/min

Injection volume: 20 μ L

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 1.5 for the hydromorphone peak

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of the labeled amount of hydromorphone hydrochloride ($C_{17}H_{19}NO_3 \cdot HCl$) in the portion of Oral Solution taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times 100$$

r_U = peak response from the *Sample solution*

r_S = peak response from the *Standard solution*

C_S = concentration of USP Hydromorphone Hydrochloride RS in the *Standard solution* (mg/mL)

C_U = nominal concentration of hydromorphone hydrochloride in the *Sample solution* (mg/mL)

Acceptance criteria: 90.0%–110.0%

IMPURITIES

ORGANIC IMPURITIES

Diluent, Solution A, and Solution B: Prepare as directed in the *Assay*.

Mobile phase: See *Table 2*.

Table 2

Time (min)	Solution A (%)	Solution B (%)
0	94	6
25	94	6
40	20	80
70	20	80
75	94	6
90	94	6

The *System suitability solution*, *Quantitation limit solution*, *Standard solution*, and *Sample solution* should be kept in a cool place, protected from light.

System suitability solution: 0.8 mg/mL of USP Hydromorphone Hydrochloride RS and 0.8 μ g/mL of USP Hydromorphone Related Compound A RS in *Diluent*

Quantitation limit solution: 0.4 μ g/mL of USP Hydromorphone Hydrochloride RS in *Diluent*

Standard solution: 4 μ g/mL of USP Hydromorphone Hydrochloride RS in *Diluent*

Sample solution: Nominally 0.4 mg/mL of hydromorphone hydrochloride in *Diluent*

Chromatographic system

(See *Chromatography* (621), *System Suitability*.)

Mode: LC

Detector: UV 220 nm

Column: 3.9-mm \times 15-cm; 5- μ m packing L1

Column temperature: 45°

Flow rate: 1 mL/min

Injection volume: 20 μ L

System suitability

Samples: *System suitability solution*, *Quantitation limit solution*, and *Standard solution*

Suitability requirements

Resolution: NLT 1.0 between the hydromorphone related compound A and hydromorphone peaks, *System suitability solution*

Signal-to-noise ratio: NLT 10:1, *Quantitation limit solution*

Tailing factor: NMT 1.5 for the hydromorphone peak, *Standard solution*

Relative standard deviation: NMT 5.0%, *Standard solution*

Analysis

Samples: *Diluent*, *Standard solution*, and *Sample solution*

Calculate the percentage of any specified or unspecified degradation product in the portion of Oral Solution taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times (1/F) \times 100$$

r_U = peak response for each degradation product found, including those in *Table 3*, from the *Sample solution*

r_S = peak response of hydromorphone from the *Standard solution*

2 Hydromorphone

C_s = concentration of USP Hydromorphone Hydrochloride RS in the *Standard solution* (mg/mL)

C_u = nominal concentration of hydromorphone hydrochloride in the *Sample solution* (mg/mL)

F = relative response factor for the corresponding individual specified or unspecified impurity (see *Table 3*)

Calculate the total degradation products by summing the percentage of all individual specified and unspecified degradation products determined to be at a level of 0.1% or greater, excluding the known process impurities, as indicated in *Table 3*.

Acceptance criteria: See *Table 3*. Disregard peaks corresponding to those from the *Diluent*, peaks that elute before a relative retention time of about 0.50, except for any peak with a relative retention time of about 0.34, and peaks that elute at the relative retention times of the process-related substances designated in *Table 3*.

Table 3

Name	Relative Retention Time	Relative Response Factor	Acceptance Criteria, NMT (%)
Specified and unidentified degradation product	0.34	1.0	0.2
8-Hydroxy-hydromorphone ^{a,b}	0.50	—	—
Dihydromorphine (DHM) ^{a,c}	0.61	—	—
Morphine ^{a,d}	0.65	—	—
Hydromorphone <i>N</i> -oxide ^{e,f}	0.79	0.87	0.2
Hydromorphone	1.0	—	—
2,2'-Bis-hydromorphone dihydrochloride ^{e,g}	2.02	1.7	0.2
Individual unspecified degradation products	—	1.0	0.2
Total degradation products	—	—	1.0

^a Process impurity.

^b 4,5 α -Epoxy-17-methylmorphinan-3,8-diol-6-one.

^c 4,5 α -Epoxy-17-methylmorphinan-3,6 α -diol.

^d 7,8-Didehydro-4,5 α -epoxy-17-methylmorphinan-3,6 α -diol.

^e Degradation product.

^f 4,5 α -Epoxy-3-hydroxy-17-methylmorphinan-6-one *N*-oxide.

^g 2,2'-Bihydromorphone.

SPECIFIC TESTS

- **MICROBIAL ENUMERATION TESTS (61) and TESTS FOR SPECIFIED MICROORGANISMS (62):** The total aerobic microbial

count does not exceed 10² cfu/mL, and the total yeasts and molds count does not exceed 10 cfu/mL. It meets the requirements of the test for the absence of *Escherichia coli*.

Change to read:

- **PH (791):** 3.5 (RB 1-Aug-2013)–6.5

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in tight containers, protected from light. Store at 25°, excursions permitted between 15° and 30°.
- **LABELING:** Identify in the product labeling any preservative used in the Oral Solution.

Change to read:

- **USP REFERENCE STANDARDS (11)**
USP Hydromorphone Hydrochloride RS
USP Hydromorphone Related Compound A RS
7,8-Didehydro-4,5 α -epoxy-3-hydroxy-17-methylmorphinan-6-one hydrochloride.
C₁₇H₁₇NO₃ 319.78 (RB 1-Aug-2013)