

[x](#)

How to Use

- **Searching:** Type keyword in search field at top of page. Search by all or part of a monograph title. For searches using multiple criteria, you will find items that match each of the specified criteria unless quotation marks are used.
 - For example, a search on Aminosalicyclic Acid Tablets will result in anything that contains “Aminosalicyclic” OR “Acid” OR “Tablets”
 - A search for “Aminosalicyclic Acid Tablets” will result in anything that specifically contains “Aminosalicyclic Acid Tablets”
- **Sorting:** Click on any column header title to sort alphabetically or chronologically in ascending or descending order. Note: the page load column is sorted alphabetically so that a number is ordered by first digit vs. by the actual number; thus, numbers will not always be in order.
 - For example, page 2178 will come before page 74 on a page sort.
- **Downloading:** You can download the Errata table in Comma-separated Value (.csv). The download will include the Errata that you have filtered on.
- **Importing:** You will need to import the file into Excel or Open Office with UTF-8 encoding, as opposed to simply opening it. To import, open Excel or Open Office and select import from the File drop-down. Depending on the version you are using, you should be presented with import formatting options to include UTF-8 as one of the first steps. Importing via UTF-8 should eliminate odd character conversions.

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
Sodium Sulfite, RE	USP36–NF31	1196	28-Mar-2014	1-Apr-2014	USP38–NF33	USP38–NF33	Line 2: Change

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
Anhydrous	AGEN TS/Reagent Specifications								[7753-83-7] to: [7757-83-7]
ONDANSETRO N ORAL SOLUTION	Assay	USP36–NF31	4586	28-Mar-2014		1-Apr-2014	USP38–NF33	USP38–NF33	Line 7 of Procedure: Change (293.36/329.82) $100(C/V)(r_U / r_S)$ to: (293.36 / 329.83) $100(C /$ $V)(r_U / r_S)$ AND Line 8 of Procedure: Change 329.82 to: 329.83
BUTYLPARAB EN	IM PUR ITIES/Related Sub stances/ Chromatographi c system	Second Supplement to USP36–NF31	6551	28-Mar-2014		1-Apr-2014	USP38–NF33	USP38–NF33	After the Column section: Add Column temperature: 35°
PHARMACEU TICAL COMPO UNDING— NONSTERILE PREPARATION	STABILITY CRITERIA AND BEYOND-USE DATIN G/General	Revision Bulletin (Official January 01, 2014)	Online	28-Mar-2014		1-Apr-2014	Second Supplement to USP37–NF32	Second Supplement to USP37–NF32	Line 7 of Paragraph 1: Change (see the General Notices

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
S	<i>Guidelines for Assigning Beyond-Use Dates</i>								<i>and Requirements, Preservation, Packaging, Storage, and Labeling)</i> to: (see <659>)
MELPHALAN TABLETS	<i>Dissolution <711></i>	USP36–NF31	4232	28-Mar-2014		1-Apr-2014	USP38–NF33	USP38–NF33	Line 1 of <i>Mobile phase</i> : Change Prepare a filtered and degassed mixture of water, acetonitrile, ammonium acetate, glacial acetic acid, and triethylamine (1 500:500:2:2:0.4). Make adjustments if necessary (see <i>System Suitability</i> under <i>Chromatography</i> <621>). to: Transfer 2 grams of

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
THIMEROSAL	IM PUR ITIES/ <i>Mercury Ions</i>	<i>USP36–NF31</i>	5368	28-Mar-2014		1-Apr-2014	<i>USP38–NF33</i>	<i>USP38–NF33</i>	<p>ammonium acetate, 2 mL of glacial acetic acid, and 0.4 mL of triethylamine to a suitable flask containing 1500 mL of water and 500 mL of acetonitrile. Stir until all solids are dissolved and well mixed, then filter and degas.</p> <p>Line 19 of <i>Analysis</i>: $C_S =$ concentration of mercuric chloride in the <i>Standard solution</i> (mg/mL) to: $C_S =$ concentration of mercuric chloride in <i>Sample solution</i></p>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
ONDANSETRO Assay N INJECTION		<i>Second Supplement to USP36–NF31</i>	Online	28-Mar-2014		1-Apr-2014	<i>USP38–NF33</i>	<i>USP38–NF33</i>	<i>B (mg/mL)</i> Line 7 of <i>Procedure:</i> Change (293.36 / 329.82)(25C / V)(<i>r</i> _U / <i>r</i> _S) to: (293.36 / 329.83)(25C / V)(<i>r</i> _U / <i>r</i> _S) AND Line 8 of <i>Procedure:</i> Change 329.82 to: 329.83
MAGNESIUM ALUMINUM SILICATE	IM PUR ITIES/ <i>Arsenic</i> , <i>Method I <211></i>	<i>USP36–NF31</i>	2073	28-Mar-2014		1-Apr-2014	<i>USP38–NF33</i>	<i>USP38–NF33</i>	Line 1 of <i>Standard preparation:</i> Change Prepare as directed in the chapter. to: Transfer 5.0 mL (5 ?g of arsenic) of the <i>Standard Arsenic Solution</i> to a 25-mL

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<p>volumetric flask, and add dilute hydrochloric acid (1:25) to volume.</p> <p>AND</p> <p>Delete:</p> <p><i>Control preparation:</i> Transfer 5.0 mL (5 ?g of As) of the <i>Standard preparation</i> to a 25-mL volumetric flask, and add dilute hydrochloric acid (1:25) to volume.</p> <p>AND</p> <p>Line 1 of <i>Acceptance criteria:</i> Change the absorbance due to any red color from the <i>Test preparation</i> does not exceed that produced by the</p>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
SODIUM ACETATE	IM PUR ITIES/ <i>Inorganic Impurities/Potassium</i>	USP36–NF31	5147	28-Mar-2014		1-Apr-2014	USP38–NF33	USP38–NF33	<p><i>Control preparation.</i></p> <p>to:</p> <p>the absorbance due to any red color from the <i>Test preparation</i> does not exceed that produced by the <i>Standard preparation.</i></p> <p>Line 1 of <i>Sample solution:</i> Change Equivalent to 600 mg/mL of anhydrous sodium acetate to:</p> <p>Dissolve the equivalent of 3 g of anhydrous sodium acetate in 5 mL of water.</p> <p>AND</p> <p>Line 1 of <i>Analysis:</i></p>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
CARBIDOPA AND LEVODOPA ORALLY DISINTEGRATING TABLETS	ASSAY	<i>Second Supplement to USP36–NF31</i>	6580	28-Mar-2014		1-Apr-2014	<i>USP38–NF33</i>	<i>USP38–NF33</i>	Change To 5 mL of <i>Sample solution</i> add to: To the <i>Sample solution</i> add Line 2 of <i>Procedure</i> : Change Inject the <i>Sample solution</i> within 2 h of preparation. Protect the volumetric solutions from light. to: Protect the volumetric solutions from light.
OIL- AND WATER-SOLUBLE VITAMINS WITH MINERALS CAPSULES	STREN GTH/ <i>Molybdenum, Method 1</i>	<i>Second Supplement to USP36–NF31</i>	6372	31-Jan-2014		1-Feb-2014	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Line 2 of <i>Standard solutions</i> : Change 2.0 to: 5.0
GLYCERYL	IDENTIFICATION	<i>USP36–NF31</i>	2033	31-Jan-2014		1-Feb-2014	<i>USP38–NF33</i>	<i>Second</i>	Line 19 of

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
TRISTEARATE	N/ <i>Fatty Acid Composition</i>							<i>Supplement to USP37–NF32</i>	<p><i>Analysis:</i> Change Result = $[(F_{MC} \times P_{FA1} \times A_{MC}) / (P_{MC} \times A_{FA1})] \times 100$ to: Result = $(F_{MC} \times P_{FA1} \times A_{MC}) / (P_{MC} \times A_{FA1})$</p>
DICLOFENAC IM SODIUM DELA PUR YED-RELEASE ITIES/ <i>Organic</i> TABLETS <i>Impurities/Procedure</i>		USP36–NF31	3221	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	<p>Line 22 of <i>Analysis:</i> Change r_S = peak response for each impurity from the <i>Standard solution</i> to: r_S = peak response of diclofenac related compound A from the <i>Standard solution</i></p>
TIZANIDINE TABLETS	IM PUR	USP36–NF31	5408	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to</i>	Line 9 of <i>Sample</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
		ITIES/ <i>Organic Impurities/Procedure</i>						USP37–NF32	<i>solution:</i> Change 45-µm or finer pore size to: 0.45-µm or finer pore size
OXALIPLATIN INJECTION	IM PURITIES/ <i>Limit of Oxalic Acid</i>	<i>First Supplement to USP36–NF31</i>	6033	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 3 of <i>Analysis:</i> Change Calculate the percentage of each impurity to: Calculate the percentage of oxalic acid
WATER-SOLUBLE VITAMINS WITH MINERALS CAPSULES	ST REN GTH/ <i>Molybdenum, Method 1</i>	<i>Second Supplement to USP36–NF31</i>	6479	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 2 of <i>Standard solutions:</i> Change 2.0 to: 5.0
CEFTAZIDIME FOR INJECTION	Assay	USP36–NF31	2887	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 6 of <i>Procedure:</i> Change 250,000[C/W (100 ? m ? s)](r _U / r _S) to:

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
RIBOFLAVIN 5 IM ?-PHOSPHATE PURITIES/ SODIUM	Free <i>Riboflavin and Riboflavin Diphosphates</i>	<i>USP36–NF31</i>	5037	31-Jan-2014		1-Feb-2014	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	25,000{C/[W (100 ? m ? s)]}(r _U / r _S) Line 1 of <i>Emission wavelength:</i> Change 530 nm to: 530 nm (monoc hromator-based detector) or 470 nm (filtered-type detector)
CYCLOMETHI CONE	ASSAY/ <i>Proce dure/System suitability</i>	<i>First Supplement to USP36–NF31</i>	5909	31-Jan-2014		1-Feb-2014	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Line 2 of <i>Suitability requirements:</i> Change <i>Relative standard deviation:</i> NMT 2.0% for cyclomethicone 4, cyclomethicone 5, and cyclomethicone 6, <i>Standard solution A, Standard solution B,</i> and

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
OIL- AND WAT ST ER-SOLUBLE REN VITAMINS GTH/ <i>Chromium</i> WITH MINERALS ORAL SOLUTION	<i>Second</i> <i>Supplement to</i> <i>USP36–NF31</i>	6399	31-Jan-2014	1-Feb-2014	<i>USP38–NF33</i>	<i>Second</i> <i>Supplement to</i> <i>USP37–NF32</i>	<i>Standard</i> <i>solution C</i> to: <i>Relative</i> <i>standard</i> <i>deviation: NMT</i> 2.0% for cyclomethicone 4, <i>Standard</i> <i>solution A</i> ; NMT 2.0% for cyclomethicone 5, <i>Standard</i> <i>solution B</i> ; NMT 2.0% for cyclomethicone 6, <i>Standard</i> <i>solution C</i> Line 13 of <i>Analysis:</i> Change C = concentration of chromium in the <i>Standard</i> <i>solution</i> (µg/mL) to: C = concentration of chromium in the <i>Sample solution</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
ANAGRELIDE HYDROCHLORIDE	ASSAY/ <i>Procedure/Chromatographic system/System suitability</i>	USP36–NF31	2500	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	(µg/mL) Line 2 of <i>Suitability requirements:</i> Change <i>Column efficiency:</i> NMT to: <i>Column efficiency:</i> NLT
DILTIAZEM HYDROCHLORIDE	IMPURITIES/ <i>Organic Impurities</i>	USP36–NF31	3258	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 21 of <i>Analysis:</i> Change r_S = response of each impurity peak from the <i>Standard solution</i> to: r_S = peak response of desacetyl diltiazem from the <i>Standard solution</i>
VANCOMYCIN HYDROCHLORIDE	SPECIFIC TESTS/ <i>Composition of Vancomycin</i>	USP36–NF31	5543	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 28 of <i>Analysis:</i> Change Result = $[(r_I / (D \times r_B) + r_A)] \times 100$

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
CONTAINERS GLASS- EVALUATION OF INNER SURFACE DURABILITY	EVALUATION OF THE INNER SURFACE D UR ABILITIES/ Aggressive Screening Conditions	Second Supplement to USP36–NF31	6221	31-Jan-2014		1-Feb-2014	USP38–NF33	Second Supplement to USP37–NF32	to: Result = $\{r_i / [(D \times r_B) + r_A]\} \times 100$ Row 1 of Column 3 of Table 4: Change 3% Citric Acid pH 8.0 to: 3% Sodium Citrate pH 8.0
ISOPROPYL ALCOHOL	ASSAY/ Procedure	Second Supplement to USP36–NF31	6638	31-Jan-2014		1-Feb-2014	USP38–NF33	Second Supplement to USP37–NF32	Line 17 of Chromatographi c system: Change Linear velocity: 35 cm/s to: Flow rate: 2.3 mL/min AND Row 2 of Column 1 of Table 2: Change Diethyl ether to: Ethyl ether
CLARITHROM	ASSAY/	USP36–NF31	3016	31-Jan-2014		1-Feb-2014	USP38–NF33	Second	Line 1 of

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
YCIN	<i>Procedure/System suitability</i>							<i>Supplement to USP37–NF32</i>	<i>Samples: Change Standard solution 2 and Standard solution 4 to: Standard solution 1, Standard solution 2, and Standard solution 4 AND Line 13 of Suitability requirements: Change Relative standard deviation: NMT 1.5%, Standard solution 2 to: Relative standard deviation: NMT 1.5%, Standard solution 1</i>
RISPERIDONE TABLETS	<i>Dissolution <711></i>	<i>USP36–NF31</i>	5065	31-Jan-2014		1-Feb-2014	<i>USP38–NF33</i>	<i>Second Supplement to</i>	<i>Line 4 of Chromatographi</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
								USP37–NF32	c system: Change Chromatograph the <i>Standard solution</i> and the <i>Test solution</i> as directed for <i>Procedure:</i> to: Chromatograph the <i>Standard solution</i> as directed for <i>Procedure:</i>
ATOMOXETIN E HYDROCHL ORIDE	IM PUR ITIES/ <i>Organic Impurities, Procedure 2</i>	<i>First Supplement to USP36–NF31</i>	5947	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 5 of <i>System suitability solution:</i> Change dissolving the Reference Standards in ethanol, to: dissolving the Reference Standards in absolute alcohol, AND Line 2 of

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
OIL- AND WATER-SOLUBLE VITAMINS WITH MINERALS ORAL SOLUTION	<i>Second Supplement to USP36–NF31</i>	6399	31-Jan-2014	1-Feb-2014	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	<p><i>Sample solution:</i> Change dissolving it in ethanol, to: dissolving it in absolute alcohol, Add the test: <i>Absence of Specified Microorganisms <2022></i>: Meet the requirements of the tests for the absence of <i>Salmonella species</i>, <i>Escherichia coli</i>, and <i>Staphylococcus aureus</i></p>
CANDESARTAN ASSAY/ NILEXETIL Procedure	<i>USP36–NF31</i>	2774	31-Jan-2014	1-Feb-2014	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	<p>Line 1 of <i>Analysis</i>: Change Titrate with 8 mL of 0.1 N to: Titrate with 0.1</p>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
FENTANYL	IM PUR ITIES/ <i>Organic Impurities/Acceptance criteria</i>	USP36–NF31	3554	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	N Footnote f of Table 2: Change N -Phenyl -N -[1-(2-phenylethyl)-4-piperidinyl]acetanilide hydrochloride, or acetyl fentanyl. to: N -(1-Phenethylpiperidin-4-yl)-N-phenylacetamide.
VANCOMYCIN SPECIFIC HYDROCHLORIDE FOR INJECTION	TESTS/ <i>Content of Vancomycin</i>	USP36–NF31	5546	31-Jan-2014		1-Feb-2014	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 27 of <i>Analysis</i> : Change Result = $[r_i / (D \times r_B) + r_A] \times 100$ to: Result = $\{r_i / [(D \times r_B) + r_A]\} \times 100$
LACTATED	<i>Definition</i>	USP36–NF31	5055	22-Nov-2013		1-Dec-2013	USP38–NF33	<i>Second</i>	Line 13:

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
RINGER'S INJECTION								<i>Supplement to USP37–NF32</i>	Change 408.0 mg of chloride to: 428.0 mg of chloride
BUFFER SOLUTIONS	<i>4. Standard Buffer Solutions/4.1 Preparation</i>	<i>Second Supplement to USP36–NF31</i>	6244	22-Nov-2013		1-Dec-2013	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Line 1 of 4. <i>Boric Acid and Potassium Chloride 0.2 M:</i> Change 12.73 g/L of boric acid to: 12.37 g/L of boric acid
DIVALPROEX SODIUM EXTENDED-RELEASE TABLETS	<i>PERFORMANCE TESTS/ Dissolution <711>/Test 5</i>	<i>Second Supplement to USP36–NF31</i>	6592	22-Nov-2013		1-Dec-2013	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Line 3 of <i>Acid stage sample solution:</i> Change suitable filter of 45-µm pore size. to: suitable filter of 0.45-µm pore size. AND Line 3 of <i>Buffer stage sample solution:</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
ALMOND OIL	SPECIFIC TESTS/ <i>Sterol Composition</i>	USP36–NF31	1877	22-Nov-2013		1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	Change suitable filter of 45- μ m pore size. to: suitable filter of 0.45- μ m pore size. Row 9 of Table 2: Change ?7-Stigmastenol ?3.0% to: ?7-Stigmasteno l ?3.0%
FLUOCINONID E TOPICAL SOLUTION	<i>Alcohol content</i>	USP36–NF31	3618	22-Nov-2013		1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	Line 1 of Standard solution: Change Dilute 20.0 mL of USP Alcohol to: Dilute 20.0 mL of alcohol
MINERAL OIL	SPECIFIC TESTS	USP36–NF31	4372	22-Nov-2013		1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	Line 1 of Viscosity—Capillary Viscometer Methods <911>: Change 34.5 and 150.0 mm ² ·s ^{?1}

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
ATROPINE SULFATE	SPECIFIC TESTS/Optical Rotation, Specific Rotation <781>	<i>First Supplement to USP36–NF31</i>	5948	22-Nov-2013		1-Dec-2013	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	to: 34.5–150.0 mm ² ·s ⁻¹ Line 1 of <i>Sample solution</i> : Change 0.1 mg/mL in water to: 0.1 g/mL of Atropine Sulfate in water
CAPSICUM OLEORESIN	SPECIFIC TESTS/Limit of Nonivamide	<i>Second Supplement to USP36–NF31</i>	6577	22-Nov-2013		1-Dec-2013	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 1 of <i>Acceptance criteria</i> : Change on the dried basis to: on the anhydrous basis
RISPERIDONE ORAL SOLUTION	ADDITIONAL REQUIREMENT S/USP Reference Standards <11>	<i>Second Supplement to USP36–NF31</i>	6690	22-Nov-2013		1-Dec-2013	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 2 of USP Risperidone Related Compounds Mixture RS: Change Contains a 98.9: 0.5: 0.30: 0.3 (area %)

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<p>mixture of four compounds: 98.9% of <i>Risperidone</i>. 0.5% of <i>Risperidone cis-N-oxide: cis</i> -3-[2-[4-(6-Fluoro-1,2-benzisoxazol-3-yl)-1-piperidinyl]ethyl]-6,7,8,9-tetrahydro-2-methyl-4H-pyrido[1,2-a]pyrimidin-4-one]. 0.3% of <i>Bicyclorisperidone: 3-(4-Fluoro-2-hydroxyphenyl)-1-[2-(6,7,8,9-tetrahydro-2-methyl-4H-pyrido[1,2-a]pyrimidin-3-yl)ethyl]-1-aza-2-az</i></p>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<p>oniabicyclo[2.2.2]oct-2-ene iodide. 0.3% of <i>Z-Oxime</i>: (<i>Z</i>)-3-[2-[4-(2,4-Difluorophenyl) (hydroxyimino)methyl]-1-piperidinylethyl]-6,7,8,9-tetrahydro-2-methyl-4<i>H</i>-pyrido[1,2-<i>a</i>]pyrimidin-4-one</p> <p>.</p> <p>to: Contains a mixture of the following four compounds: 98.9% of <i>Risperidone</i>. 0.5% of <i>Risperidone cis-N-oxide</i>: <i>cis</i>-3-[2-[4-(6-Fluoro-1,2-benzisoxazol-3-yl)-1-piperidinylethyl]-6,7,</p>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<p>8,9-tetrahydro-2-methyl-4H-pyrido[1,2-a]pyrimidin-4-one-N-oxide.</p> <p>0.3% of <i>Bicyclorisperidone</i>: 3-(4-Fluoro-2-hydroxyphenyl)-1-[2-(6,7,8,9-tetrahydro-2-methyl-4-oxo-4H-pyrido[1,2-a]pyrimidin-3-yl)ethyl]-1-aza-2-azoniabicyclo[2.2.2]oct-2-ene iodide.</p> <p>0.3% of <i>Z-Oxime</i>: (Z)-3-[2-[4-(2,4-Difluorophenyl) (hydroxyimino)methyl]-1-piperidiny</p>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
BOSWELLIA SERRATA	COMPOSITION <i>/Content of Keto-Derivatives of ?-Boswellic Acids/System suitability/Suitability requirements</i>	USP36–NF31	1366	22-Nov-2013		1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	Line 1 of <i>Tailing factor</i> . Change 11-keto-?-boswellic acid peak to: 3-acetyl-11-keto-?-boswellic acid peak Change the subsection head <i>Capillary zone electrophoresis system</i> (see <i>Capillary Electrophoresis under Biotechnology-Derived Articles—Test <1047></i>)— to: <i>Capillary zone</i>
APROTININ	<i>Limit of des-Ala-aprotinin and des-Ala-des-Gly-aprotinin</i>	USP36–NF31	2522	22-Nov-2013		1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	Change the subsection head <i>Capillary zone electrophoresis system</i> (see <i>Capillary Electrophoresis under Biotechnology-Derived Articles—Test <1047></i>)— to: <i>Capillary zone</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
KETOROLAC T IM ROMETHAMIN PUR E TABLETS ITIES/ <i>Organic Impurities</i>	USP36–NF31	4042	22-Nov-2013	1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	<i>electrophoresis system—</i> Line 1 of <i>Sample solution:</i> Change Proceed as directed for the <i>Sample stock solution</i> in the Assay. to: Proceed as directed for the <i>Sample solution</i> in the Assay.
RISPERIDONE <i>USP Reference standards <11></i>	USP36–NF31	5063	22-Nov-2013	1-Dec-2013	USP38–NF33	Second Supplement to USP37–NF32	Line 7 of USP Risperidone System Suitability Mixture RS: Change 9-Hydroxyrisperidone-(<i>6RS</i>)-3-[2-[4-(6-fluoro-1,2-benzisoxazol-3-yl)piperidin-1-yl]ethyl]-2,6-dimethyl-6,7,8,9-tet

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
									rahydr o-4H -pyri do[1,2-a]pyrimidin-4-one . to: 9-Hydroxyrispe ridone: (9RS)-3-{2-[4-(6-fluor o-1,2-benzisoxa zol-3-yl)piperidi n-1-yl]ethyl}-9-h ydroxy-2-methyl -6,7,8,9-tetrahy dro-4H -pyri do[1,2-a]pyrimidin-4-one .
CHINESE SALVIA	COMPOSITION	<i>Second Supplement to /Content of Salvianolic Acid USP36–NF31 B</i>	6331	22-Nov-2013		1-Dec-2013	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Line 16 of <i>Analysis:</i> Change <i>W</i> = weight of Chinese Salvia used to prepare the <i>Sample solution</i> (mg) to: <i>W</i> = weight of

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
LORATADINE	IM PURITIES/Organic Impurities, Procedure 1	<i>Second Supplement to USP36–NF31</i>	6650	22-Nov-2013		1-Dec-2013	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Chinese Salvia used to prepare the <i>Sample stock solution</i> (mg) Line 3 of Note: Change 4,8-dichloro-6,11-dihydro-5H-benzo[5,6]cyclohepta[1,2-b]pyridin-11-one to: 4,8-dichloro-5,6-dihydro-11H-benzo[5,6]cyclohepta[1,2-b]pyridin-11-one
POLYETHYLENE GLYCOL MONOMETHYL ETHER	IM PURITIES/Limit of 2-Methoxyethanol	<i>USP36–NF31</i>	2142	22-Nov-2013		1-Dec-2013	<i>USP38–NF33</i>	<i>Second Supplement to USP37–NF32</i>	Line 11 of <i>Calibration</i> : Change On the two <i>Calibration</i> plots, to: On the <i>Calibration</i> plot,

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
GADOTERIDOL INJECTION	<i>Bacterial endotoxins <85></i>	USP36–NF31 3701	22-Nov-2013	1-Dec-2013	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 2: Change gadoteridol. to: Gadoteridol Injection.
OCTINOXATE	<i>Assay</i>	USP36–NF31 4556	22-Nov-2013	1-Dec-2013	USP38–NF33	<i>Second Supplement to USP37–NF32</i>	Line 3 of <i>Chromatographic system:</i> Change a 0.32-mm x 25-m column that contains coating G1, to: a 0.32-mm x 25-m column with 0.25-µm thickness of phase G1 coating,

Pagination

- [First page « First](#)
- [Previous page ‹ Previous](#)
- ...
- [Page 30](#)
- [Page 31](#)
- [Page 32](#)
- [Page 33](#)
- [Page 34](#)

-
- [Page 35](#)
 - [Page 36](#)
 - [Page 37](#)
 - [Page 38](#)
 - ...
 - [Next page Next ›](#)
 - [Last page Last »](#)