



Metoprolol Succinate Extended-Release Tablets

Type of Posting	Notice of Intent to Revise
Posting Date	30-Aug-2024
Targeted Official Date	To Be Determined, Revision Bulletin
Expert Committee	Small Molecules 2

In accordance with the Rules and Procedures of the Council of Experts and the [Pending Monograph Guideline](#), this is to provide notice that the Small Molecules 2 Expert Committee intends to revise the Metoprolol Succinate Extended-Release Tablets monograph.

Based on the supporting data received from a manufacturer awaiting FDA approval, the Expert Committee proposes to add *Dissolution Test 9*. The revision also necessitates a change in the table numbering in the test for *Organic Impurities*.

- The analytical procedure in *Dissolution Test 9* was validated using the Agilent ZORBAX SB C8 brand of column with L7 packing. The typical retention time for metoprolol succinate is about 2.5 min.

The proposed revision is contingent on FDA approval of a product that meets the proposed monograph specifications. The proposed revision will be published as a Revision Bulletin and an official date will be assigned to coincide as closely as possible with the FDA approval of the associated product.

See below for additional information about the proposed text.¹

Should you have any questions, please contact Brice Wagner, Scientist III (301-998-6832 or brice.wagner@usp.org).

¹ This text is not the official version of a *USP–NF* monograph and may not reflect the full and accurate contents of the currently official monograph. Please refer to the current edition of the *USP–NF* for official text.

USP provides this text to indicate changes that we anticipate will be made official once the product subject to this proposed revision under the Pending Monograph Program receives FDA approval. Once FDA approval is granted for the associated revision request, a Revision Bulletin will be posted that will include the changes indicated herein, as well as any changes indicated in the product's final approval, combined with the text of the monograph as effective on the date of approval. Any revisions made to a monograph under the Pending Monograph Program that are posted without prior publication for comment in the *Pharmacopeial Forum* must also meet the requirements outlined in the [USP Guideline on Use of Accelerated Processes for Revisions to the USP–NF](#).

Metoprolol Succinate Extended-Release Tablets

DEFINITION

Metoprolol Succinate Extended-Release Tablets contain NLT 90.0% and NMT 110.0% of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$.

IDENTIFICATION

- **A.** The retention time of the major peak of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in the *Assay*.
- **B.** The UV spectrum of the major peak of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in the *Assay*.

ASSAY

• PROCEDURE

Buffer: Mix 50 mL of 1 M [monobasic sodium phosphate](#) and 8.0 mL of 1 M [phosphoric acid](#), and dilute with [water](#) to 1000 mL. If necessary, adjust with 1 M [monobasic potassium phosphate](#) or 1 M [phosphoric acid](#) to a pH of 3.0.

Mobile phase: [Acetonitrile](#) and *Buffer* (25:75)

Standard solution: 0.05 mg/mL of [USP Metoprolol Succinate RS](#) in *Mobile phase*

Sample stock solution: Nominally 1 mg/mL of metoprolol succinate prepared as follows. Transfer a suitable number of Tablets to a suitable volumetric flask, add about 5 mL of [water](#), and allow the Tablets to disintegrate. Add a volume of [alcohol](#) to fill 30% of the flask volume, and shake for 30 min. Add a portion of 0.1 N [hydrochloric acid](#) to fill 50% of the flask volume, and shake for an additional 30 min. Dilute with 0.1 N [hydrochloric acid](#) to volume. Filter, and discard the first 10 mL of the filtrate.

Sample solution: Nominally 0.05 mg/mL of metoprolol succinate from the *Sample stock solution* in *Mobile phase*

Chromatographic system

(See [Chromatography](#) (621), [System Suitability](#).)

Mode: LC

Detector: UV 280 nm. For *Identification B*, use a diode array detector in the range of 200–400 nm.

Column: 4-mm × 12.5-cm; 5- μ m packing [L7](#)

Flow rate: 1 mL/min

Injection volume: 40 μ L

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2.0

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ in the portion of Tablets taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times 100$$

r_U = peak response of metoprolol from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* (mg/mL)

C_U = nominal concentration of metoprolol succinate in the *Sample solution* (mg/mL)

Acceptance criteria: 90.0%–110.0%

PERFORMANCE TESTS

Change to read:

- **DISSOLUTION** <711>

Test 1

Medium: pH 6.8 phosphate buffer (see [Reagents, Indicators, and Solutions—Buffer Solutions](#)); 500 mL

Apparatus 2: 50 rpm

Times: 1, 4, 8, and 20 h

Buffer and Mobile phase: Prepare as directed in the Assay.

Standard solution: A known concentration of [USP Metoprolol Succinate RS](#) in *Medium*

Sample solution: Pass a portion of the solution under test through a suitable filter.

Chromatographic system

(See [Chromatography](#) <621>, [System Suitability](#).)

Mode: LC

Detector: UV 280 nm

Column: 4-mm × 12.5-cm; 5- μ m packing [L7](#)

Flow rate: 1 mL/min

Injection volume: 40 μ L

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2.0

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ dissolved at each time point.

Tolerances: See [Table 1](#).

Table 1

Time (h)	Amount Dissolved (%)
1	NMT 25
4	20–40
8	40–60
20	NLT 80

The percentages of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ dissolved at the times specified conform to [Dissolution](#) (711), [Acceptance Table 2](#).

Test 2: If the product complies with this test, the labeling indicates that the product meets USP *Dissolution Test 2*.

Medium: Simulated gastric fluid without enzyme, pH 1.2; 500 mL

Apparatus 2: 75 rpm

Times: 1, 4, 8, and 20 h

Buffer: 1 M [monobasic sodium phosphate](#), 1 M [phosphoric acid](#), and [water](#) (50:8:942). If necessary, adjust with 1 M [monobasic sodium phosphate](#) or 1 M [phosphoric acid](#) to a pH of 3.0.

Mobile phase: [Acetonitrile](#) and *Buffer* (250:750)

Standard solution: Prepare a solution of [USP Metoprolol Succinate RS](#) in *Medium* as directed in [Table 2](#).

Table 2

Tablet Strength (mg)	Concentration (mg/mL)
200	0.380
100	0.190
50	0.095
25	0.048

Sample solution: Pass the solution under test through a suitable filter.

Chromatographic system

(See [Chromatography](#) (621), [System Suitability](#).)

Mode: LC

Detector: UV 280 nm

Column: 4.0-mm × 12.5-cm; 4-μm packing [L7](#)

Flow rate: 1 mL/min

Injection volume: See [Table 3](#).

Table 3

Tablet Strength (mg)	Volume (µL)
25	40
50	20
100	10
200	5

Run time: NLT 2 times the retention time of metoprolol

System suitability

Sample: *Standard solution*

Suitability requirements

Column efficiency: NLT 1500 theoretical plates

Tailing factor: NMT 2.0

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the concentration (C_i) of metoprolol succinate dissolved in *Medium* at each time point

(i):

$$\text{Result} = (r_U/r_S) \times C_S$$

r_U = peak response of metoprolol from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* (mg/mL)

Calculate the percentage of the labeled amount of metoprolol succinate [$(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4$] dissolved (Q_i), at each time point (i):

$$\text{Result}_1 = C_1 \times V \times (1/L) \times 100$$

$$\text{Result}_2 = \{[C_2 \times (V - V_S)] + (C_1 \times V_S)\} \times (1/L) \times 100$$

$$\text{Result}_3 = (\{C_3 \times [V - (2 \times V_S)]\} + [(C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

$$\text{Result}_4 = (\{C_4 \times [V - (3 \times V_S)]\} + [(C_3 + C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

C_i = concentration of metoprolol succinate in the portion of sample withdrawn at time point (i) (mg/mL)

V = volume of *Medium*, 500 mL

L = label claim (mg/Tablet)

V_S = volume of the *Sample solution* withdrawn from the *Medium* (mL)

Tolerances: See [Table 4](#).

Table 4

Time Point (i)	Time (h)	Amount Dissolved (%)
1	1	NMT 20
2	4	20–40
3	8	55–85
4	20	NLT 80

The percentages of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ dissolved at the times specified conform to [Dissolution](#) (711), [Acceptance Table 2](#).

Test 4: If the product complies with this test, the labeling indicates that the product meets USP *Dissolution Test 4*.

Medium: Phosphate buffer, pH 6.8 (dissolve 6.8 g of [monobasic potassium phosphate](#) and 0.93 g of [sodium hydroxide](#) in 1 L of [water](#); adjust with a [sodium hydroxide](#) solution to a pH of 6.8); 500 mL

Apparatus 2: 50 rpm

Times: 1, 4, 8, and 24 h

Buffer: 5.0 mL/L of [triethylamine](#) in [water](#). Adjust with [phosphoric acid](#) to a pH of 3.0.

Mobile phase: [Methanol](#) and *Buffer* (40:60)

Standard solution: Prepare a solution of [USP Metoprolol Succinate RS](#) in *Medium* as directed in [Table 5](#).

Table 5

Tablet Strength (mg)	Concentration (mg/mL)
200	0.4
100	0.2
50	0.1
25	0.05

Sample solution: Withdraw a 10-mL aliquot at each time point. Pass the solution under test through a suitable filter of 0.45- μ m pore size. Replace the portion withdrawn with an equal volume of *Medium*.

Chromatographic system

(See [Chromatography](#) (621), [System Suitability](#).)

Mode: LC

Detector: UV 223 nm

Column: 4.6-mm \times 25-cm; 5- μ m packing [L1](#)

Column temperature: 30°

Flow rate: 1.5 mL/min

Injection volume: 5 µL

Run time: NLT 2 times the retention time of metoprolol

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2.0

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the concentration (C_i) of metoprolol succinate dissolved in *Medium* at each time point (i):

$$\text{Result} = (r_U/r_S) \times C_S$$

r_U = peak response of metoprolol from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* (mg/mL)

Calculate the percentage of the labeled amount of metoprolol succinate [$(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4$] dissolved (Q_i), at each time point (i):

$$\text{Result}_1 = C_1 \times V \times (1/L) \times 100$$

$$\text{Result}_2 = [(C_2 \times V) + (C_1 \times V_S)] \times (1/L) \times 100$$

$$\text{Result}_3 = \{(C_3 \times V) + [(C_2 + C_1) \times V_S]\} \times (1/L) \times 100$$

$$\text{Result}_4 = \{(C_4 \times V) + [(C_3 + C_2 + C_1) \times V_S]\} \times (1/L) \times 100$$

C_i = concentration of metoprolol succinate in the portion of sample withdrawn at time point (i) (mg/mL)

V = volume of *Medium*, 500 mL

L = label claim (mg/Tablet)

V_S = volume of the *Sample solution* withdrawn from the *Medium* (mL)

Tolerances: See [Table 6](#).

Table 6

Time Point (i)	Time (h)	Amount Dissolved (Tablets labeled 25 mg) (%)	Amount Dissolved (Tablets labeled 50, 100, and 200 mg) (%)
1	1	NMT 20	NMT 20
2	4	20–40	15–35
3	8	42–67	38–64
4	24	NLT 80	NLT 80

The percentages of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ dissolved at the times specified conform to [Dissolution <711>](#), [Acceptance Table 2](#).

Test 5: If the product complies with this test, the labeling indicates that the product meets USP *Dissolution Test 5*.

Medium: Phosphate buffer, pH 6.8 (dissolve 27.22 g of [monobasic potassium phosphate](#) and 3.6 g of [sodium hydroxide](#) in 4 L of water; adjust with 1 N [sodium hydroxide](#) or [phosphoric acid](#) to a pH of 6.8); 500 mL

Apparatus 2: 50 rpm, with sinkers

Times: 1, 4, 8, and 20 h

Buffer: Transfer 3.0 mL of [triethylamine](#) and 1.0 mL of [phosphoric acid](#) to a 1000-mL volumetric flask that contains 600 mL of [water](#). Dilute with [water](#) to volume.

Mobile phase: [Acetonitrile](#) and *Buffer* (25:75)

Standard solution: Prepare a solution of [USP Metoprolol Succinate RS](#) in *Medium* as directed in [Table Z](#).

Table 7

Tablet Strength (mg)	Concentration (mg/mL)
200	0.2
100	0.2
50	0.05
25	0.05

Sample solution: Withdraw a 10-mL aliquot at each time point. Pass the solution under test through a suitable filter of 0.45- μ m pore size. Replace the portion withdrawn with an equal volume of *Medium*.

Chromatographic system

(See [Chromatography <621>](#), [System Suitability](#).)

Mode: LC

Detector: UV 280 nm

Column: 4.6-mm \times 15-cm; 5- μ m packing [L7](#)

Column temperature: 40°

Flow rate: 1.5 mL/min

Injection volume: 40 μ L for 25 and 50 mg; 10 μ L for 100 and 200 mg

Run time: NLT 2 times the retention time of metoprolol

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2.0

Relative standard deviation: NMT 3.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the concentration (C_i) of metoprolol succinate dissolved in *Medium* at each time point

(i):

$$\text{Result} = (r_U/r_S) \times C_S$$

r_U = peak response of metoprolol from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* (mg/mL)

Calculate the percentage of the labeled amount of metoprolol succinate [$(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4$] dissolved (Q_i), at each time point (i):

$$\text{Result}_1 = C_1 \times V \times (1/L) \times 100$$

$$\text{Result}_2 = [(C_2 \times V) + (C_1 \times V_S)] \times (1/L) \times 100$$

$$\text{Result}_3 = \{(C_3 \times V) + [(C_2 + C_1) \times V_S]\} \times (1/L) \times 100$$

$$\text{Result}_4 = \{(C_4 \times V) + [(C_3 + C_2 + C_1) \times V_S]\} \times (1/L) \times 100$$

C_i = concentration of metoprolol succinate in the portion of sample withdrawn at time point (i) (mg/mL)

V = volume of *Medium*, 500 mL

L = label claim (mg/Tablet)

V_S = volume of the *Sample solution* withdrawn from the *Medium* (mL)

Tolerances: See [Table 8](#).

Table 8

Time Point (i)	Time (h)	Amount Dissolved (%)
1	1	NMT 10
2	4	5–30
3	8	30–55
4	20	NLT 75

The percentages of the labeled amount of metoprolol succinate [$(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4$] dissolved at the times specified conform to [Dissolution](#) (711), [Acceptance Table 2](#).

Test 7: If the product complies with this test, the labeling indicates that the product meets USP *Dissolution Test 7*.

Medium: Phosphate buffer, pH 6.8 (dissolve 6.8 g of [monobasic potassium phosphate](#) and 0.9 g of [sodium hydroxide](#) in 1 L of [water](#); adjust with 1 N [sodium hydroxide](#) or 1 M [phosphoric acid](#) to a pH of 6.8), deaerated; 500 mL

Apparatus 2: 50 rpm

Times: 1, 4, 8, 12, and 24 h for Tablets labeled 25 and 50 mg and 1, 4, 8, and 20 h for Tablets labeled 100 and 200 mg

Buffer: Mix 50 mL of 1 M [monobasic sodium phosphate dihydrate](#) and 8.0 mL of 1 M [phosphoric acid](#), and dilute with [water](#) to 1000 mL. Adjust with 1 M [monobasic sodium phosphate dihydrate](#) or 1 M [phosphoric acid](#) to a pH of 3.0.

Mobile phase: [Acetonitrile](#) and *Buffer* (22:78)

Standard solution: Known concentrations of [USP Metoprolol Succinate RS](#) in *Medium* are listed in [Table 9](#).

Table 9

Tablet Strength (mg)	Concentration (mg/mL)
100	0.2
50/200	0.1
25	0.05

Sample solution: Withdraw an 8-mL aliquot of the solution under test at each time point. For Tablets labeled 25, 50, and 100 mg, no further dilution is required. For Tablets labeled 200 mg, transfer 5.0 mL of the aliquot withdrawn to a 20-mL volumetric flask, and dilute with *Medium* to volume. Centrifuge, if needed. Pass the supernatant through a suitable filter of 0.45- μ m pore size.

Chromatographic system

(See [Chromatography](#) (621), [System Suitability](#).)

Mode: LC

Detector: UV 280 nm

Column: 4.6-mm \times 15-cm; 5- μ m packing [L7](#)

Column temperature: 30°

Flow rate: 0.9 mL/min

Injection volume: 10 μ L

Run time: NLT 2.6 times the retention time of metoprolol

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2.0

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the concentration (C_i) of metoprolol succinate dissolved at each time point (i):

$$\text{Result} = (r_U/r_S) \times C_S \times D$$

r_U = peak response of metoprolol from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* (mg/mL)

D = dilution factor for the *Sample solution*

Calculate the percentage of the labeled amount of metoprolol succinate [$(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4$] dissolved (Q_i), at each time point (i):

$$\text{Result}_1 = C_1 \times V \times (1/L) \times 100$$

$$\text{Result}_2 = \{[C_2 \times (V - V_S)] + (C_1 \times V_S)\} \times (1/L) \times 100$$

$$\text{Result}_3 = (\{C_3 \times [V - (2 \times V_S)]\} + [(C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

$$\text{Result}_4 = (\{C_4 \times [V - (3 \times V_S)]\} + [(C_3 + C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

$$\text{Result}_5 = (\{C_5 \times [V - (4 \times V_S)]\} + [(C_4 + C_3 + C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

C_i = concentration of metoprolol succinate in the portion of sample withdrawn at time point (i) (mg/mL)

V = volume of *Medium*, 500 mL

L = label claim of metoprolol succinate (mg/Tablet)

V_S = volume of the sample withdrawn at each time point, 8 mL

Tolerances: See [Table 10](#) and [Table 11](#).

Table 10

Time Point (i)	Time (h)	Amount Dissolved (Tablets labeled 25 and 50 mg) (%)
1	1	NMT 15
2	4	10-30
3	8	35-55
4	12	55-75
5	24	NLT 80

Table 11

Time Point (i)	Time (h)	Amount Dissolved (Tablets labeled 100 and 200 mg) (%)
1	1	NMT 17
2	4	17-37

Time Point (i)	Time (h)	Amount Dissolved (Tablets labeled 100 and 200 mg) (%)
3	8	42–62
4	20	NLT 80

The percentages of the labeled amount of metoprolol succinate [$(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4$] dissolved at the times specified conform to [Dissolution <711>](#), [Acceptance Table 2](#).

▲Test 9: If the product complies with this test, the labeling indicates that the product meets USP [Dissolution Test 9](#).

Medium: Phosphate buffer, pH 6.8; 500 mL

Apparatus 2: 50 rpm

Times: 1, 4, 8, 20, and 24 h

Buffer: Mix 50 mL of 1 M [monobasic sodium phosphate](#) and 8.0 mL of 1 M [phosphoric acid](#), and dilute with [water](#) to 1000 mL. Adjust with 1 M [monobasic sodium phosphate](#) or 1 M [phosphoric acid](#) to a pH of 3.0.

Mobile phase: [Acetonitrile](#) and [Buffer](#) (25:75)

Standard solution: 0.05 mg/mL of [USP Metoprolol Succinate RS](#) in [Medium](#)

Sample solution: At the times specified, withdraw a portion of the solution under test, and pass the solution under test through a suitable filter of 0.45- μ m pore size.

Chromatographic system

(See [Chromatography <621>](#), [System Suitability](#).)

Mode: LC

Detector: UV 280 nm

Column: 4-mm \times 12.5-cm; 5- μ m packing [L7](#)

Flow rate: 1 mL/min

Injection volume: 10 μ L

Run time: NLT 2 times the retention time of metoprolol

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2.0

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the concentration (C_i) of metoprolol succinate dissolved at each time point (i):

$$\text{Result} = (r_U/r_S) \times C_S$$

r_U = peak response of metoprolol from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* (mg/mL)

Calculate the percentage of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ dissolved at each time point (*i*):

$$\text{Result}_1 = C_1 \times V \times (1/L) \times 100$$

$$\text{Result}_2 = \{[C_2 \times (V - V_S)] + (C_1 \times V_S)\} \times (1/L) \times 100$$

$$\text{Result}_3 = (\{C_3 \times [V - (2 \times V_S)]\} + [(C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

$$\text{Result}_4 = (\{C_4 \times [V - (3 \times V_S)]\} + [(C_3 + C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

$$\text{Result}_5 = (\{C_5 \times [V - (4 \times V_S)]\} + [(C_4 + C_3 + C_2 + C_1) \times V_S]) \times (1/L) \times 100$$

C_i = concentration of metoprolol succinate in the portion of sample withdrawn at time point (*i*) (mg/mL)

V = volume of Medium, 500 mL

L = label claim of metoprolol succinate (mg/Tablet)

V_S = volume of the sample withdrawn at each time point (*i*) (mL)

Tolerances: See [Table 12](#).

Table 12

Time Point (<i>i</i>)	Time (h)	Amount Dissolved (%)
1	1	NMT 20
2	4	20-40
3	8	40-60
4	20	68-88
5	24	NLT 75

The percentages of the labeled amount of metoprolol succinate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_4]$ dissolved at the times specified conform to [Dissolution <711>](#), [Acceptance Table 2](#). ▲ (TBD)

- **UNIFORMITY OF DOSAGE UNITS** <905>: Meet the requirements

IMPURITIES

Change to read:

- **ORGANIC IMPURITIES**

Buffer: 1.15 mL of [phosphoric acid](#) in 2 L of [water](#). Add 2.6 g of [sodium dodecyl sulfate](#). Sonicate to dissolve.

Solution A: [Methanol](#) and *Buffer* (30:70)

Solution B: [Acetonitrile](#) and *Buffer* (75:25)

Mobile phase: See [Table 13](#).

Table 13 (TBD)

Time (min)	Solution A (%)	Solution B (%)
0	65	35
20	65	35
25	40	60
30	35	65
35	35	65
37	65	35
50	65	35

Diluent: [Acetonitrile](#) and *Buffer* (40:60)

System suitability solution: 3 µg/mL of [USP Metoprolol Related Compound A RS](#) and 1 mg/mL of [USP Metoprolol Succinate RS](#) in *Diluent*

Standard solution: 3 µg/mL of [USP Metoprolol Succinate RS](#) in *Diluent*

Sensitivity solution: 0.5 µg/mL of [USP Metoprolol Succinate RS](#) from *Standard solution* in *Diluent*

Sample solution: Nominally 1 mg/mL of metoprolol succinate from Tablets prepared as follows. Transfer a portion of finely powdered Tablets (NLT 20), equivalent to 50 mg of metoprolol succinate, to a 50-mL volumetric flask. Add *Diluent* to fill 60% of the flask volume and sonicate for 30 min with intermittent shaking. Dilute with *Diluent* to volume. Pass the solution through a suitable filter of 0.45-µm pore size.

Chromatographic system

(See [Chromatography](#) (621), [System Suitability](#).)

Mode: LC

Detector: UV 223 nm

Column: 4.6-mm × 15-cm; 5-µm packing [L1](#)

Column temperature: 30°

Flow rate: 1 mL/min

Injection volume: 10 µL

System suitability

Samples: *System suitability solution*, *Standard solution*, and *Sensitivity solution*

[NOTE—The relative retention times for succinic acid (the counter ion), metoprolol related compound A, and metoprolol are about 0.1, 0.83, and 1.0, respectively.]

Suitability requirements

Resolution: NLT 2.0 between metoprolol related compound A and metoprolol, *System suitability solution*

Relative standard deviation: NMT 5.0%, *Standard solution*

Signal-to-noise ratio: NLT 10, *Sensitivity solution*

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of any unspecified impurity in the portion of Tablets taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times 100$$

r_U = peak response of any unspecified impurity from the *Sample solution*

r_S = peak response of metoprolol from the *Standard solution*

C_S = concentration of [USP Metoprolol Succinate RS](#) in the *Standard solution* ($\mu\text{g/mL}$)

C_U = nominal concentration of metoprolol succinate in the *Sample solution* ($\mu\text{g/mL}$)

Acceptance criteria: The reporting threshold is 0.05%.

Any unspecified impurity: NMT 0.20%

Total impurities: NMT 0.75%

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in tight containers, and store at controlled room temperature.
- **LABELING:** Label it to indicate the content of metoprolol succinate and its equivalent, expressed as metoprolol tartrate $[(C_{15}H_{25}NO_3)_2 \cdot C_4H_6O_6]$. When more than one *Dissolution* test is given, the labeling states the *Dissolution* test used only if *Test 1* is not used.
- **USP REFERENCE STANDARDS** [\(11\)](#)

[USP Metoprolol Related Compound A RS](#)

1-Ethylamino-3-[4-(2-methoxyethyl)phenoxy]propan-2-ol.

$C_{14}H_{23}NO_3$ 253.34

[USP Metoprolol Succinate RS](#)

Page Information:

Not Applicable

Current DocID:

© The United States Pharmacopeial Convention *All Rights Reserved.*