













- [WATER DETERMINATION \(921\), Method I](#)

**Sample:** 2.0 g

**Acceptance criteria:** NMT 1.0%

#### **ADDITIONAL REQUIREMENTS**

- **PACKAGING AND STORAGE:** Preserve in tight containers, protected from direct sunlight. Avoid exposure to excessive heat.
- **LABELING:** Label it to indicate its polydispersity ( $M_w/M_n$ ) or its polydispersity range. Label it to indicate the name and amount of any added antioxidant.
- **USP REFERENCE STANDARDS (11).**
  - [USP Diethylene Glycol RS](#)
  - [USP Ethylene Glycol RS](#)
  - [USP Polyethylene Glycol 3350 RS](#)

**Polyethylene Glycol 3350** —see [Polyethylene Glycol 3350 General Monographs](#)

---

<sup>1</sup> Available from Sigma-Aldrich Corporation, or equivalent.

<sup>2</sup> Available from Sigma-Aldrich Corporation, or equivalent.

<sup>3</sup> Millipore® Millex® LCR HPLC syringe filters with hydrophilic PTFE membrane is suitable, or any other equivalent filter.

#### **Page Information:**

Not Applicable

#### **DocID:**

© 2020 The United States Pharmacopeial Convention *All Rights Reserved.*